# amoena

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**Thank you for choosing a quality product from Amoena Medizin-Orthopädie-Technik GmbH. For more than 45 years, Amoena has been synonymous with innovative and high-quality products that support the needs of women following breast surgery.**

#### **Intended Purpose /Indication**

Target group: Women after various breast surgeries

Compression textiles are designed to keep the operated breast or breasts immobile and/or stabilised during the healing phase.

Compression textiles

- can support wound healing after breast surgery.
- can stabilise the shape and function of the operated breast.
- can protect sensitive scar tissue and reduce irritation.
- can minimise the risk of postoperative oedema.
- can reduce postoperative swelling/pain.

#### **Contraindications/Risks**

Relative contraindications: In the case of the following clinical symptoms, it is up to your doctor to decide whether a compression textile can be used:

- Skin diseases
- Very sensitive scar area in the surgical zone
- Radiodermatitis
- Known material allergy (contact dermatitis)

Side effects: If the wrong size is selected, this may lead to localised pressure hypersensitivity and disorders of blood flow and/or lymphatic flow. In rare cases allergic reactions may occur. A contraindication exists if there is a known allergy to one of the components that may come into contact with the skin.

#### **General Care Instructions**

- The care instructions can be found on the care label of the compression textile.
- Before washing, close the zip and VELCRO® fastenings as well as hooks and eyes.
- It is recommended that the textile is washed in a laundry bag.
- If the compression textile is washed with other textiles, it is recommended that similar colours should be washed together.

#### **SEAMLESS COMPRESSION TEXTILES**

Art. No.	Article	Art. No.	Article
44604	Leyla ZIP-ST-LOW (rose nude)	45008	Sina HE-ST-MEDIUM (black)
44605	Leyla ZIP-ST-LOW (white)	45009	Pamela ZIP-ST-HIGH (rose nude)
44606	Leyla ZIP-ST-LOW (black)	45010	Pamela ZIP-ST-HIGH (white)
45006	Sina HE-ST-MEDIUM (rose nude)	45011	Pamela ZIP-ST-HIGH (black)
45007	Sina HE-ST-MEDIUM (white)		

#### **Instructions for Use of Seamless Compression Textiles**

To ensure an optimal healing process, sufficient immobilisation and stabilisation of the treated breast tissue is highly advisable. Optimal fixation/stabilization of the operated breast tissue is also guaranteed with double sizes and universal cups. The compression textile in the primary packaging is prepared for direct clinical use in the operating room and can be adjusted to fit the patient directly after surgery by trained personnel. Depending on the individual patient's medical indication, the Seamless Compression Bra can be combined with any compression belt.

#### **Determining the Correct Bra Size**

The patient's underbust circumference when standing, is determined by using a standard proprietary metric measuring tape. Using the table below, the underbust measurement obtained, can be allocated to an underbust size. The unique seamless knitting technology, can provide universal cup sizes and individually adjusts to the breast size and shape, from A – F cup.

Underbust measurement	68 – 77	78–87	88–97	98–107	108–117	118–127
Size	70/75	80/85	90/95	100/105	110/115	120/125

**Material Composition of Seamless Compression Textiles**

Leyla ZIP-ST-LOW	Polyamide, cotton, elastane
Sina HE-ST-MEDIUM	Polyamide, elastane
Pamela ZIP-ST-HIGH	Polyamide, elastane

For percentages of material composition, please read the care label.  
Free of latex, silicone and nickel.

**Special Notes Seamless Compression Textiles**

Do not use if packaging is damaged

**COMPRESSION TEXTILES**

Art. No.	Article
0776N	Compression Belt (white)
0778	Sarah SB (white)
0779	Sarah SB (black)
0789N	Compression Belt (black)
45043	Compression Belt (rose nude)
2161	Theraport Post Surgery Bra (white)
45044	Anatomical Belt (white)
45045	Anatomical Belt (black)
45046	Anatomical Belt (rose nude)

**Instructions for Use of Compression Textiles****Determining the Correct Bra Size**

The dimensions of the underbust and overbust should be determined using a tape measure with the patient in a standing position. The underbust measurement gives the width of the underbust while the overbust measurement gives the cup size. Depending on the medical indication, the compression bra can be combined with any compression belt.

**Theraport Post Surgery Bra instructions for use****Putting on the bra:**

Place the bra around your body with the Velcro closure facing forward. Fasten the Velcro closure underneath your chest. When it is properly in place, roll the bra up over your chest. Slip your arms into the straps and pull them up over your shoulders.

**To prepare the bra for brachytherapy:**

Find the catheter outlet location in the firmer bra material under the arm and mark it on the inside using a felt-tip pen. Then cut a small opening the size of the catheter into the material. Make sure that it is above the underbust band. Once the catheter has been fitted, you can adjust the bra and guide the end of the catheter out through the opening in the bra.

**Material Composition of Compression Textiles**

Sarah SB	Cotton, polyamide, elastane
Compression Belt	Polyamide, elastane, polyester, cotton
Theraport	Polyamide, cotton, elastane

For percentage of material composition see care label  
Free of latex, silicone and nickel.

### Storage advice and durability

A maximum shelf life of 36 months applies. The storability of compression textiles is 6 months and depends on the correct handling such as following care instructions, donning and doffing (or "dressing and undressing the compression product") and the usage of the products to a normal extent. The useful life of the medical device is shown by the hourglass symbol on the packaging. When properly treated, the compression textiles retain their medicinal effect over the full period of 6 months. If your body dimensions should change - due to your medical analysis and individual life circumstances, a new measurement and a new fitting is advisable even before the 6-month period has expired

### Special Instructions

Protect from sunlight

Store in a dry place

Single patient – multiple uses

### Disposal Instructions

Dispose of the products in household waste.

### Reporting Obligation

If you have any complaints about the product, please contact your medical retailer directly. It is only necessary to notify the manufacturer and the competent authority of the member state in the event of any serious incidents that could lead to a significant deterioration of health or to death. Serious incidents are defined in Article 2 No. 65 of Regulation (EU) 2017/745 (MDR).

### Warranty

The statutory warranty obligation applies.

### Service Information

Do you have any questions?

If so, we would like to help – please call us or write to us.  
You can find our address on page 95.

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Leyla ZIP-ST-LOW	Poliammide, cotone, elasthan
Sina HE-ST-MEDIUM	Poliammide, elasthan
Pamela ZIP-ST-HIGH	Poliammide, elasthan

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Art. n.	Articolo
0776N	Compression Belt (white)
0778	Sarah SB (white)
0779	Sarah SB (black)
0789N	Compression Belt (black)
45043	Compression Belt (rose nude)
2161	Reggiseno post-operatorio Theraport 2161 (white)
45044	Anatomical Belt (white)
45045	Anatomical Belt (black)
45046	Anatomical Belt (rose nude)

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Amoena is a registered trademark of Amoena Medizin-Orthopädie-Technik GmbH.
© 2022 All rights reserved Amoena Medizin-Orthopädie-Technik GmbH,
Amoena Sweden AB
Amoena UK Ltd