amoena

RECOVERY CARE **CURASUPPORT COMPRESSION**

Gebrauchsanweisung. Instructions for use. Mode d'emploi.
Instrucciones de uso. Instruções de uso. Istruzioni per l'uso.
Gebruiksaanwijzing. Brugsanvisning. Instruksjoner for bruk.
Användningsinstruktioner. Instrukcja użycia. Käyttöohjeet.
Návod k použití. Инструкции по применению.
Οδηγίες χρήσης. 使用説明. 使用説明書. مداخلت سالا تامي لعنت

Deutsch	DE	4
English	ENG	8
Français	FR	12
Español	ES	17
Italiano	IT	22
Nederlands	NL	26
Dansk	DK	30
Norsk	NO	34
Svenska	SE	38
Polski	PL	42
Português	PT	46
Suomalainen	FI	50
Český	CZ	54
Русский	RU	58
Ελληνικά	GR	62
中國	CN	66
日本の	JP	70
74	ARAB	العربية

POST-OP BRAS

Art. No.	Article
44604	Leyla ZIP-ST-LOW (Rose Nude)
44605	Leyla ZIP-ST-LOW (White)
44606	Leyla ZIP-ST-LOW (Black)
45006	Sina HE-ST-MEDIUM (Rose Nude)
45007	Sina HE-ST-MEDIUM (White)
45008	Sina HE-ST-MEDIUM (Black)
45009	Pamela ZIP-ST-HIGH (Rose Nude)
45010	Pamela ZIP-ST-HIGH (White)
45011	Pamela ZIP-ST-HIGH (Black)

POST-OP BRAS WITH HIGH COTTON CONTENT

Art. No.	Article
0778	Sarah SB (White)
0779	Sarah SB (Black)
2161	Theraport (White)

COMPRESSION BELTS

Art. No.	Article
0776N	Compression Belt ¹ (White)
0789N	Compression Belt ¹ (Black)
45043	Compression Belt ¹ (Rose Nude)
45044	Anatomical Compression Belt ² (White)
45045	Anatomical Compression Belt ² (Black)
45046	Anatomical Compression Belt ² (Rose Nude)

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Leyla ZIP-ST-LOW: Polyamid, Baumwolle, Elastan

Sina HE-ST-MEDIUM: Polyamid, Elastan

Pamela ZIP-ST-HIGH: Polyamid, Elastan

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Thank you for choosing a quality product from Amoena Medizin-Orthopädie-Technik GmbH. Since 1975, Amoena has stood for innovative and high-quality products for the special needs of women who have undergone breast surgery.

Intended Purpose /Indication

Target group: Women after various breast surgeries

Compression textiles are designed to keep the operated breast or breasts immobile and/or stabilised during the healing phase.

Compression textiles

- can support wound healing after breast surgery.
- can stabilise the shape and function of the operated breast.
- can protect sensitive scar tissue and reduce irritation.
- can minimise the risk of postoperative oedema.
- can reduce postoperative swelling/pain.

Contraindications/Risks

Relative contraindications: In the case of the following clinical symptoms, it is up to your doctor to decide whether a compression textile can be used:

- Skin diseases
- Very sensitive scar area in the surgical zone
- Radiodermatitis
- Known material allergy (contact dermatitis)

Side effects

If the wrong size is selected, this may lead to localised pressure hypersensitivity and disorders of blood flow and/or lymphatic flow. A contraindication exists if there is a known allergy to one of the components that may come into contact with the skin. In rare cases, residual risks such as allergic reactions may occur. Please talk to your GP if you experience any significant side effects.

General Care Instructions

- The care instructions can be found on the care label of the compression textile.
- Before washing, close the zip and VELCRO® fastenings as well as hooks and eyes.
- It is recommended that the textile is washed in a laundry bag.
- If the compression textile is washed with other textiles, it is recommended that similar colours should be washed together.

POST-OP BRA

Material Composition

Leyla ZIP-ST-LOW: Polyamide, cotton, elastane

Sina HE-ST-MEDIUM: Polyamide, elastane

Pamela ZIP-ST-HIGH: Polyamide, elastane

For percentages of material composition, please read the care label.

Free of latex, silicone and nickel.

Instructions for Use of Seamless Compression Textiles

To ensure an optimal healing process, sufficient immobilisation and stabilisation of the treated breast tissue is highly advisable. Optimal fixation/stabilization of the operated breast tissue is also guaranteed with double sizes and universal cups. The compression textile in the primary packaging is prepared for direct clinical use in the operating room and can be adjusted to fit the patient directly after surgery by trained personnel. Depending on the individual patient's medical indication, the Seamless Compression Bra can be combined with any compression belt.

Determining the Correct Bra Size

The patient's underbust circumference when standing, is determined by using a standard proprietary metric measuring tape. Using the table below, the underbust measurement obtained, can be allocated to an underbust size. The unique seamless knitting technology, can provide universal cup sizes and individually adjusts to the breast size and shape, from A – F cup.

Underbust measurement	58–67	68–77	78–87	88–97	98–107	108–117	118–127	128–137
Size	60/65	70/75	80/85	90/95	100/105	110/115	120/125	130/135

Special Notes Seamless Compression Textiles

Do not use if packaging is damaged

POST-OP BRAS WITH HIGH COTTON CONTENT / COMPRESSION BELT

Material Composition

Sarah SB: Cotton, polyamide, elastane

Anatomical Compression Belt/Compression Belt: Polyamide, elastane, polyester, cotton

Theraport: Polyamide, cotton, elastane

For percentage of material composition see care label

Free of latex, silicone and nickel.

Instructions for Use of post-op bras with high cotton content

Determining the Correct Bra Size

The dimensions of the underbust and overbust should be determined using a tape measure with the patient in a standing position. The underbust measurement gives the width of the underbust while the overbust measurement gives the cup size. Depending on the medical indication, the compression bra can be combined with any compression belt.

Theraport Post Surgery Bra instructions for use

Putting on the bra:

Place the bra around your body with the Velcro closure facing forward. Fasten the Velcro closure underneath your chest. When it is properly in place, roll the bra up over your chest. Slip your arms into the straps and pull them up over your shoulders.

To prepare the bra for brachytherapy:

Find the catheter outlet location in the firmer bra material under the arm and mark it on the inside using a felt-tip pen. Then cut a small opening the size of the catheter into the material. Make sure that it is above the underbust band. Once the catheter has been fitted, you can adjust the bra and guide the end of the catheter out through the opening in the bra.

Instructions for use of compression belt

The compression belts are fastened with the help of the hook and loop strap above the chest:

1. Before putting it on, make sure that the Amoena woven label should be readable from the front.
2. Secure the middle of the belt with the hook and loop strap, starting at the bottom of the hook and loop area in the back of the compression bra.
3. Bring the two ends of the belt around to the front and underneath the armpits.
4. Close the hook and loop strap above the chest.
5. Make sure that the compression belt exerts a comfortable and clearly noticeable pressure. If you experience pain or any other discomfort, loosen the hook and loop strap

Storage advice and durability

A maximum shelf life of 36 months and a duration of use of 6 months apply. The maximum life span is therefore limited to 42 months. The storability of compression textiles is 6 months and depends on the correct handling such as following care instructions, donning and doffing (or "dressing and undressing the compression product") and the usage of the products to a normal extent. The useful life of the medical device is shown by the hourglass symbol on the packaging. When properly treated, the compression textiles retain their medicinal effect over the full period of 6 months. If your body dimensions should change - due to your medical analysis and individual life circumstances, a new measurement and a new fitting is advisable even before the 6-month period has expired

Special Instructions

Page 78

Disposal Instructions

Dispose of the products in household waste.

Reporting Obligation

If you have any complaints about the product, please contact your medical retailer directly. It is only necessary to notify the manufacturer and the competent authority of the member state in the event of any serious incidents that could lead to a significant deterioration of health or to death. Serious incidents are defined in Article 2 No. 65 of Regulation (EU) 2017/745 (MDR).

Warranty

The statutory warranty obligation applies.

Service Information

Do you have any questions?

If so, we would like to help – please call us or write to us.

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Amoena is a registered trademark of Amoena Medizin-Orthopädie-Technik GmbH.
© 2023 All rights reserved Amoena Medizin-Orthopädie-Technik GmbH,
Amoena Sweden AB
Amoena UK Ltd