

amoena

Hygienically prepared for direct clinical use.

# LEYLA

## POST-SURGICAL COMPRESSION BRA

Designed to support and optimize surgical results for various indications:

- Breast Reconstruction
- Breast Conserving Surgery/Lumpectomy
- Mastectomy
- Breast Reduction
- Breast Augmentation
- Breast Lifting
- Oncoplastic Breast Surgeries



**Hygienically prepared for direct clinical use** through a prewashing and disinfection process

**Sufficient immobilization and stabilization** of the surgical site and breast tissue encourages healing

**Minimized discomfort on surgical site** through seamless, gentle fabric that reduces skin irritation and unnecessary pressure

**Easy fitting** due to seamless technology and color-coded sizing



- 1 **Seamless design** using strong yet gentle fabric provides low compression and feels soft and comfortable on sensitive skin
- 2 **Targeted stabilization zones** assist with optimal lymphatic flow around the torso
- 3 **Adjustable wide lined straps** with hook and loop fasteners for easy fitting and wearing comfort
- 4 **Zip front closure** for easy access to the surgical site with an interior hook and eye placket for added security
- 5 **Soft fabric flap** at bottom of zipper protects sensitive skin from chaffing
- 6 **Seamless cup pockets** can hold a post-surgical form securely in place
- 7 **Belt accessory** (sold separately) additionally stabilizes the breast area, attaching to the hook and loop fastener flap in the back of the garment to help keep it in place
- 8 **Also available in Black**, to help hide stains

**Style** 44605  White  
44606  Black

**Material** 49% Nylon, 43% Cotton, 8% Spandex  
*Free of latex, silicone and nickel.*

**Insurance** PDAC registration pending



**Color-coded sizing helps you easily identify the right size from your stock.**

Underbust Measurement					
<b>Centimeters</b>	68-77	78-87	88-97	98-107	108-117
<b>Inches</b>	27-30	31-34	35-38	39-42	43-46
	↓	↓	↓	↓	↓
<b>Size</b>	32/34	36/38	40/42	44/46	48/50
<b>Cup Size</b>	Universal (One Size)				
<b>Symbol</b>	○	■	△	●	□

Take a moment to find out more and watch an informational video at [amoena.com/us-en/Leyla](http://amoena.com/us-en/Leyla)