

amoena



NEW

LEYLA
HYGIENICALLY
PREPARED
POST-SURGICAL
COMPRESSION BRA

PREWASHED & PACKAGED FOR DIRECT CLINICAL USE

DESIGNED TO SUPPORT AND OPTIMIZE SURGICAL RESULTS FOR VARIOUS INDICATIONS




- Breast Reconstruction
- Breast Conserving Surgery/Lumpectomy
- Mastectomy
- Breast Reduction
- Breast Augmentation
- Breast Lifting
- Oncoplastic Breast Surgeries

LEYLA POST-SURGICAL COMPRESSION BRA

Style	44605	<input type="checkbox"/> White
	44606	<input type="checkbox"/> Black
Sizes	32/34; 36/38; 40/42; 44/46; 48/50	
Cup Size	Universal (One Size)	
Material	49% Nylon, 43% Cotton, 8% Spandex <i>Free of latex, silicone and nickel.</i>	
Insurance	PDAC registration pending	

Color-coded sizing helps you easily identify the right size from your stock. [See back cover for details.](#)

FEATURES



Prewashed and packaged for direct clinical use in the OR

Designed to securely hold a Compression Belt in place

Also available in Black, to help hide stains



CHOOSE LEYLA FOR A SUPERIOR OUTCOME



- **Hygienically prepared for direct clinical use** through a prewashing and disinfection process
- **Sufficient immobilization and stabilization** of the surgical site and breast tissue encourages healing
- **Minimized discomfort on surgical site** through seamless, gentle fabric that reduces skin irritation and unnecessary pressure
- **Easy fitting** due to seamless technology and color-coded sizing

HOW TO FIND THE RIGHT SIZE

Color-coded sizing helps you easily identify the right size from your stock.

- 1 While wearing a bra, have your patient stand up with her arms at her side.
- 2 Measure her **underbust** by measuring around her torso where the bra band sits. Make sure the measuring tape is placed firmly and securely against her body.
- 3 Find the **underbust measurement** in the chart below to find the size and corresponding color-coded symbol.

UNDERBUST MEASUREMENT					
Centimeters	68-77	78-87	88-97	98-107	108-117
Inches	27-30	31-34	35-38	39-42	43-46
Size	32/34	36/38	40/42	44/46	48/50
Symbol	○	■	△	●	□



Take a moment to watch an informational video at amoena.com/us-en/Leyla

WE'RE HERE FOR YOU.

CONTACT YOUR AMOENA ACCOUNT MANAGER FOR MORE INFORMATION.

CUSTOMER SERVICE

T 1-800-741-0078

E USAmoenaCSR@amoena.com

ONLINE ORDERING

b2b.amoena.com/us-en

Amoena USA Corporation
1701 Barrett Lakes Blvd., Ste. 410, Kennesaw, GA USA 30144
amoena.com

amoena
Supporting Confidence